

Phone: 2538-3683 (Direct)

2538-2664 / 0203 : 2538-6442

Fax: 2538-6442 E-mail: madhyamgrammunicipality@vsnl.net

# **MADHYAMGRAM MUNICIPALITY**

P.O.: Madhyamgram
Dist.: North 24 Parganas, Kolkata - 700129

## NOTICEINVITING-eQuotation

## NIQNO:WBMAD/MM/NIQ-03e/2024-25

## **NIQNO.**MM/EO/1340/2024-25

Dated:06/08/2024

On behalf of The Executive Officer, Madhyamgram Municipality, invites e Quotations (electronic tender process) from agencies/consultant having knowledge for quoting the following materials are invited for submission of quotation. **Rates quoted should be inclusive of all taxes.** 

Sl.No.	Descriptionofsupplyof materials	Unit	Qty	Rate
1	Quotation for preparationof estimates of Supply and delivery at site of a Diathermy Machine for Gyno / General OT for Madhyamgram Municipality.	No.	1	
	Specification:			
	MONOPOLARCUT:PURESINEWAVE390KHz			
	LOW- CONTINUES, 300 W AT 300 OHMS CF 1.5,LOW VOLT			
	PURE-CONTINUES,400WAT300OHMMS,CF1.5,HIGHVOLT BLEND			
	- BURSTS 29 KHz, 250 W AT 3000HMS CF 2.5			
	ENDOCUT-PULSE,400WAT300OHMS,CF1.5			
	BIPOLARCUT:PURESINEWAVE390KHz			
	MACRO- CONTINUES OUTPUT, 100W AT 100 W AT 1000HMS CF 1.5, HIGH VOLT, RAPID			
	PULSEBICUT-PULSEOUTPUT,100WAT100OHMSCF1.5,HIGHVOLT			
	LEAKAGECURRENTS			
	.RFLEAKAGECURRENT.MONOPOLAR:LESSTHAN150maRMS			
	.BIPOLAR:LESSTHAN60MaRMS			
	.Lowfrequencyleakagecurrents			
	.NORMALPOLARITY,INTACTGROUND:LESSTHAN10Ua			
	.NORMALPOLARITY,groundopen:lessthan50Ua			
	FEEDBACKSYSTEM			
	.3Microprocessorcontrolled6sensetechnologywithsamplingrateabove4000sense/second			
	.Maximumoutputpowerdeviationlessthan15%or5watts,whicheverisgreateratany setting with per-			
	98			
	INPUTPOWER			
	Operatingvoltage:180-265voltAC,50-60Hz Max. power input :1000 VA, fuse : 6A, UPS Back up :1.2KVA			
	MONOPOLARCOAG: DAMPED SINE WAVE 460 KHz			
	SOFT-CPNTINUESAT84KHz,120WAT500OHMS,CF6.2			
	FULGURATE-REPETITIONFREQUENCY30KHz,120WAT500OHMS,CF7.0			
	SPRAY- RANDOMIZED REPETITION 30 KHz LESS THAN f LESS THAN 50KHz, 120W			
	AT 500 OHMS CF 8.0			
	BIPOLARCOAG:PURESINEWAVE390KHz			
	MICRO- CONTINUES OUTPUT, 100W AT 100 OHMS, CF 1.5, LOW VOLT			
	STANDERD-CONTINUESOUTPUT,100WAT100OHMS,CF1.5,MEDIUMVOLT			
	FORCE- CONTINUES OUTPUT, 100W AT 100 OHMS, CF 1.5, HIGH VOLT <b>SEALER</b> : PURE SINE WAVE 390KHz,			
	Rapidseal-pulseoutputmode,3level,100Wat100ohms,CF1.5 PREMSYSTEM			
	MEASURINGFREQUENCY:80KHz(+,-)10KHz MEASURING CURRENT: LESS THAN 10 uA			
	Acceptableresistancerange			
	DualareaPREMpatientreturnelectrode:5–135ohms(adaptive)			
	SINGLEAREAPATIENTRETURNELECTRODE:LESSTHAN10OHMS			
	DESIGNSTANDARD			
	CLASSIFICATIONACCORDINGTOTHEEC-DIRECTIVE93/42/EEC:CLASSIIb GENERAL			
	SAFETY STANDERDS : IEC 60601-1-1, IEC 60601-2-2			
	EMI/EMCCOMPATIBILITYSTANDERD:IEC60601-1-2 IS			
	SAFETY STANDARD : IS 7583			
	DIMENSIONS			
	LXWXHINCM-38X31X12,WEIGHINKGLESSTHAN6			

Sl.No.	Descriptionofsupplyof materials	Unit	Qty	Rate
1		No.	1	
	FEATURESOFDIATHERMY	110.		
	5INCHLCDTOUCHSCREEN:ALLOWSQUICK,EASY&INTERCTIVE EXPERIAENCE TO			
	USERS.			
	6SENSE TECHNOLOGY : MAINTAINS SET POWER OVER WIDER RANGE OF			
	TISSUE FEEDBACK TECHNOLOGY.			
	DDEM DATIENTEDETUDN ELECTRODEMONITODINGDROVIDEGGAFETY ACAIGT			
	PREM: PATIENTRETURN ELECTRODEMONITORINGPROVIDESSAFETY AGAIST THE PATIENT RETURN ELECTRODE SITE BURNS.			
	VESSEL SEALING MODE: SEAL UPTO 7 MM TISSUE BUNDLES & VESSELS WITH			
	REUSABLE VESSEL ESEALING ATTACHMENTS.			
	AUTO STOPS: AUTOMATICALLY RF CURRENT AFTER COMPLETION OF VESSEL			
	SEALIN.			
	ENDO CUT MODE: ALLOWS PULSE INTERVAL CONTROLLED CUTTING FOR			
	POLYPECTOMY & PAPILLOTOMY.			
	RANDOMIZED SPRAY COAG : PROVIDES HOMOGENIOUS, SUPERFICIAL &			
	FASTER COAGULATION OVER LARGER AREA.			
	PULSE BICUT: FOR MODERN BIPOLARCUTTING APPLICATION WITH MINIMAL			
	THERMAL DAMAGE.  BIPOLAR CUT & COAG MODES: INDEPENDENT DISPLAY SETTING USING ITS			
	RESPECTIVE FOOTSWITCH.			
	BICOAG ALARM: AUDIO FEEDBACK AFTER COMPLETION OF BIPOLAR COAG			
	REDUCES CHARRING & STICKING.			
	FOOTSWITCH: ENHANCED WITH A TOGGLE BUTTON TO SWITCH BETWEEN			
	MODES ALONG WITH ON SCREEN INDICATION.UNIT ALSO SUPPORTS WIRELESS			
	FOOTSWITCH FOR CLEAN WIRE FREE SURGICAL ENVIRONMENT.			
	SMART SURGICAL ASSIST. : THE PROGRAM MANAGER ALLOWS SMART			
	ACCESS TO PROGRAMS & SETTINGS.			
	SELF- ILLUMINATING SOCKET: HELPS DURING POWER SETTINGS &RF			
	ACTIVATION TO INDICATE USE OF RESPECTIVE TO INDICATE USE OF RESPECTIVE RF PORT.			
	SELF DIAGONOSIS: ENSURES THE ESU WITH ACCESSORIES ARE SAFE FOR			
	FUNCTIONING & DISPLAY ERROR CODES IF ANY.			
	USER PROGRAMMABLE: 30 USER SETTABLE PROGRAMS WITH PROCEDURE &			
	SURGEON NAME.			
	REMOTE POWER SETTINGS ( OPTIONAL) : TO CHANGE MODE AND POWER			
	FROM STERILE FIELD BY NORMAL HANDSWITCH PENCIL.			
	TIME OUT FUNCTION ( OPTIONAL): AUTOMATICALLY STOPS THE HF OUTPUT			
	INCASE OF PROLONGED ACTIVATION.  LIDCRADARH ITV . ARGON REAM CONGULATION LILTRASONOC ASPIRATOR			
	UPGRADABILITY: ARGON BEAM CONGULATION, ULTRASONOC ASPIRATOR, SMOKE EVACUTOR, SUCTIONIRRGATION MODULE.			
	SURGICAL APPLICATION: SUITABLE FOR ALL KIND OF OPEN, ENDOSCOPIC			
	LAPROSCOPIC & MINIMALLY INVASIVE SURGERYSUCH AS PLASTIC, NEURO,			
	ENT, SPINE , OPTHAL, ORTHO, GENERAL SURGERY, PEDIATRIC, ONCO,			
	GYNOCOLOGY, HEMIA GL ENDOCUT, PAPILLOTOMY, POLYPECTOMY, UNDER			
	WATER, UROLOGY GYNEC, TCRE, HYSTEROSCOPY, THORACTOMY, THORASIC&			
	CARDIAC SURGERY.			
	VESSEL SEALING APPLICATIONS: LAPROSCOPIC/ RADICAL / VAGINAL /			
	ABDOMINAL HYSTERECTOMY, APPENDECTOMY, COLECTOMY, CYSTECTOMY,			
	GASTTRIC BYPASS, HEMORHOIDECMY, LAVH , LIVER RESECTION, NEDHDECTOMY NISSEN FUNDODLICATION SDLENECTOMY			
	NEPHRECTOMY, NISSEN FUNDOPLICATION, SPLENECTOMY.			

Table-2:-	
Dateofuploadingofe-NIQ,and Quotation Documents online (PublishingDate)	07/08/2024at 11.00AM
Documentdownloaded/salestart date(on line)	07/08/2024at 11.00AM
Quotationsubmissionclosingdate(online)	17/08/2024upto4.00PM
Date&timeofopeningtender(Technicalbid)	20/08/2024at4.00PM

	The tender shall remain open for acceptance till 90 days from the date of opening bid and accepted rate shall remain valid		
Validityof offer			
	duringtheservicesexecutionperiodincludingextended period		
	ifany.		
	1. Experience Certificate		
Eli aihility Cuitania	2. ValidTradeLicence		
EligibilityCriteria	3. PanCard&ValidGst Certificate		
	4. LatestP.TaxCert.		
	5.ITReturn(AY-2023-2024or2024-2025)		

- ➤ Intending Quotationers/ Suppliers can search and download NIe-Q and other Quotation documents electronically by logging on to the website <a href="https://wbtenders.gov.in">https://wbtenders.gov.in</a>, using his Digital Signature Certificate (DSC). This is the mode of collection of Quotation documents electronically.
- ${\color{blue} \blacktriangleright} \ \ For participating in the Quotation, the Quotationer do not pay any amount as an Earnest Money.$

### Notes:-

- A. Submissionofbidbythebidder,doesnotconstituteany guarantee for opening of pricebids ince opening of "PriceBid" depends of fulfilment of eligibility criteria for the tender which shall be decided by the Accepting Officer after verification of documents submitted by the bidder along with "Technical Bid" to be opened on the time and date as mentioned above.
- B. TheentirequotationtobecarriedoutasperinstructionofTIA.

Executive Officer

Madhyamgram Municipality
Executive Officer

Madhyamgram Municipality
North 24 Parganas